

Copyright

Our company owns all rights of this unpublished work and intends to maintain it as a confidential work. We may also seek to maintain this work as an unpublished copyright. This publication is to be used solely for the purpose of reference or operation of our software system. No part of this work can be disseminated for other purposes.

In the event of inadvertent or deliberate publication, we intend to enforce its right to this work under copyright laws as a published work. Those having access to this work may not copy, use, or disclose the information in this work unless expressly authorized by our company.

All information contained in this publication is believed to be correct. We shall not be liable for errors contained herein nor for incidental or consequential damages in connection with the furnishing, performance, or use of this material. The information this publication refers to is protected by copyrights or patents and does not convey any license under the patent rights of our company, nor the rights of others. We do not assume any liability arising out of any infringements of patents or other rights of third parties.

Content of the manual is subject to change without prior notice.

Version: V1.0I-P

Issue Date: August 12, 2010

CONTENT

CHAPTER 1 Introduction	1
1.1 Brief Introduction	1
1.2 Safety Information	1
1.3 Electromagnetic Compatibility	3
1.4 Equipment Classification	3
1.5 Accessories	3
1.6 Equipment symbol	4
CHAPTER 2 BASIC OPERATION	5
2.1 Outer View	5
2.2 Install the batteries	5
2.3 Connect the sensor	6
2.4 SpO ₂ Monitoring	6
2.5 Factors that may affect the measurement	7
CHAPTER 3 DETAILED OPERATION	9
3.1 Display	9
3.2 Key Functions	10
3.3 Mode Introduction	10
3.4 Data replay and transmission	11
3.5 Alarm	12
CHAPTER 4 MAINTENANCE AND REPAIR	13
4.1 Maintenance	13
4.2 Troubleshooting	13
4.3 Warranty and Repair	14
APPENDIX A Specifications	17
APPENDIX B The declare of manufacturer	19

CHAPTER 1 INTRODUCTION

1.1 Brief Introduction

Thank you for purchasing the handheld pulse oximeter for the functional oxygen saturation of arterial hemoglobin (SpO₂) and Pulse Rate (PR) measurement. The handheld pulse oximeter features PR tone modulation, data storage and data transmission capabilities. Please read the operator's manual carefully before using this instrument.

Notice: The illustrations used in this manual may differ slightly from the appearance of the actual product.

Intended use: The pulse oximeter is a portable non-invasive device that intended for spot checking, displaying, storing and transmitting the oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate of single adult and pediatric patient in hospital (including clinical use in surgeries, anesthesia, intensive care and etc.) . Not for continuous monitoring.


Measurement principle

The principle of pulse oximetry is based on the red and infrared (IR) light absorption of oxygenated and deoxygenated hemoglobin present in the circulating blood. Oxygenated hemoglobin absorbs more IR and allows more red light to pass through. Deoxygenated hemoglobin conversely absorbs more red light and allows IR light to pass through. The detector probe is placed on the finger. The probe contains two light emitting diodes (LED's), one in the visible red spectrum (660nm) and one in the IR spectrum (940 nm). The beams of light from this probe pass through the tissues and some light is absorbed by the blood and soft tissues depending on hemoglobin concentration. The amount of light absorption at each light frequency is dependent on the degree of oxygenation of hemoglobin within the tissues.

The microprocessor can select out the absorbance of the pulsatile fraction of blood, i.e. that due to arterial blood, from constant absorbance due to non-pulsatile venous or capillary blood and other tissue pigments.

Principle of Operation

1. Red and Infrared-ray Emitter Diode
2. Red and Infrared-ray Receptor Diode



1.2 Safety Information

Warning, Precaution and Notice

Warning, Precaution and Notice in the manual are special information that prompts the operator's attention.

Warning - Information concerning something that could possibly hurt the patient or operator.

Precaution - Reminds the user to pay close attention to device operation, failure of which may cause abnormal function of the instrument.

Notice - Informs the user of other important information by suggestion, requirement and supplement.

Warnings

- Please read this manual carefully before using this device. The user must check that the equipment functions safely and ensure that it is in proper working condition before being used.
- Do not use the handheld pulse oximeter in an explosive atmosphere.
- Do not use the handheld pulse oximeter in an MRI or CT environment.
- The handheld pulse oximeter has no SpO₂ alarms and PR alarms; it is not for continuous monitoring, as indicated by the symbol.
- The handheld pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
- Check the handheld pulse oximeter sensor application site frequently to determine the positioning of the sensor and circulation and skin sensitivity of the patient.
- Do not stretch the adhesive tape while applying the pulse oximeter sensor. This may cause inaccurate readings or skin blisters.
- Prolonged use of the probe/sensor or the patient's condition may require changing the sensor site periodically. Change the sensor site and check skin integrity, circulatory status, and correct alignment at least every 4 hours. Prolonged use may cause blisters, skin deterioration, and discomfort.
- Sensor malfunction may cause inaccurate data possibly resulting in patient injury or death, so pay close attention to the sensor and inspect it often.
- Worn-out data cables may also cause inaccurate data, so if the data is used as a reference to treat a patient, pay special attention to data cable and check it more frequently.
- Do not tangle the SpO₂ cable with the wires of ES (Electrosurgery) equipment.
- Single use accessories should never be reused.

Precautions

- Operation of the handheld pulse oximeter may be affected by the use of an electrosurgical unit (ESU).
- The handheld pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO₂ measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO₂ measurement.
- Do not sterilize the device using autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid which may cause inaccurate readings. The device is not intended for sterilization.
- Unplug the sensor from the monitor before cleaning or disinfecting it.
- If liquid is accidentally spilled on the unit, clean and dry it thoroughly before reuse.
- Do not try to apply the SpO₂ and NIBP measurement on the same arm at the same time. This could potentially affect measurement accuracy.

Notices

- Operation of this device in an electromagnetic field may influence its accuracy.
- SpO₂ measurements may be influenced by high ambient light, especially sunlight. Shield the sensor area if necessary.
- Dyes introduced into the bloodstream, such as methylene blue, indocyanine green, indigo carmine, and fluorescein, may influence the accuracy of the SpO₂ readings.
- Any condition that restricts blood flow, such as use of a blood pressure cuff or extremes in systemic vascular resistance, may cause a failure to determine accurate pulse rate and SpO₂ readings.

- Remove fingernail polish or artificial fingernails before applying SpO₂ sensors. Fingernail polish or artificial fingernails may cause inaccurate SpO₂ readings.
- Optical cross talk can occur when two or more sensors are located in adjoining areas. It can be eliminated by covering each site with opaque material. Optical cross talk may adversely affect the accuracy of the SpO₂ readings.
- Obstructions or dirt on the sensor's red light or detector may cause a sensor failure. Make sure there are no obstructions and the sensor is clean.
- For routine equipment maintenance, please refer to the service procedures at the associated section as indicated in the manual.
- As to the other concerns for attention, please carefully look through the specific chapter in this instruction.

1.3 Electromagnetic Compatibility

This oximeter is designed and tested in compliance with the EMC standard, complying with the international standard for the EMC of the medical electrical device - IEC 60601-1-2. However, because of the proliferation of radio frequency transmitting equipment and other sources of electrical noise in the health-care and home environments (e.g. cellular phones, mobile two-way radios, electrical appliances) it is possible that high levels of such interference due to close proximity or strength of a source, may result in disruption of performance of this device.

This apparatus complies with the IEC 60601-1-2 international standard. The requirements of this international standard are: CISPR11, GROP1, and CLASS B.

1.4 Equipment Classification

Classification according to IEC-60601-1	
According to the type of protection against Electrical shock:	Internal electrical power source equipment
According to the degree of protection against Electrical shock:	Type BF equipment
According to the degree of protection against harmful ingress of water.	IPX1
According to the methods of sterilization or disinfection	Non-sterilizable: Use of Liquid surface disinfectants only.
According to the mode of operation:	Continuous operation/ Spot-checking
The equipment is not suitable for use in the presence of a flammable anesthetic mixture air or with oxygen or nitrous oxide.	

1.5 Accessories

Standard accessories:

- Operator's manual
- Fingertip sensor for adult:
Model: M-50E13CSO
- 2×1.5V AAA-Size Alkaline batteries (nanfu)
- MedView software CD
- Data cable

Optional accessories:

- Binding sensor for pediatric (15-45Kg)
Model: M-50C

Confirm that the items listed are packed with the handheld pulse oximeter. If any item on this list is missing or damaged, contact your distributor. Contact the carrier immediately if the shipping carton is damaged.

1.6 Equipment symbol

Symbols	Definition
	Attention! Refer to the relevant prompt. Read the operator's manual carefully before using the oximeter.
	Type BF applied part
	Date of Manufacture
	Low power indicator
IPX1	Protected against dripping water
	European union approval
SN	Serial number
	Not for continuous monitoring
	Prevent from rain
	Storage temperature and relative humidity

CHAPTER 2 BASIC OPERATION

2.1 Outer View




Fig.2.1 Front and rear panel

Description of Fig.2.1:

1. The interface for connecting the sensor and data transmission
2. Displaying screen
3. Power/function key
4. Setting key
5. Clip for fixing the device
6. The cover of battery compartment

2.2 Install the batteries

The oximeter can be powered by 2 AAA-Size alkaline batteries which will typically provide 16 hours of continuous operation,(by represent of nanfu for example).

When battery power is lower than $2.3 \pm 0.1V$, the sign  will flicker in its display area. Replace the battery as soon as possible. If the battery power is lower than $2.2 \pm 0.1V$, the device will power off automatically.

Pull the cover as the indication by the arrowhead. And install the batteries as the polarity markings(+ and -). Refer to fig.2.2

Be sure not to insert the batteries in wrong polarities, otherwise, the damage may be caused to the device.




Fig.2.2

2.3 Connect the sensor

Connect the oximeter sensor to sensor interface at the top of the oximeter. Ensure that the sensor is firmly plugged in.

The interface is also used for transferring data to the software for management. To transfer data, please connect the data cable to the interface. For detailed information, refer to MedView software instruction manual.

Notes:

1. Please make the raised lump side of sensor cable upside when inserting the sensor.
2. Do not insert the sensor forcefully.
3. Because the probe socket and USB socket share a common interface so you cannot take a measurement when transmitting data.

2.4 SpO₂ Monitoring

Clip the sensor to the patient's finger, and ensure that the patient's nail surface is facing upward, as shown in Fig.2.3.




Fig.2.3 Placement of the finger

2.4.1 Power on the oximeter

Press the function key (left key) to power it on. Several seconds later, the measurement value will appear.

Notice: To maintain the highest degree of accuracy, it is recommended that the finger and the oximeter sensor/probe is kept as still as possible.

2.4.2 Brightness adjustment

When you press the function key (left key) for more than one second, the brightness level headed with "Br" will be shown on the top right of the screen. You can adjust the brightness by degrees by pressing the setting key (right key). There are 7 levels of brightness. The default is level three.




Fig. 2.4

2.4.3 Mode switch

After turning on the oximeter, each time the function key (left key) is pressed, the oximeter will switch to another display mode, as shown in Fig.3.1.

2.5 Factors that may affect the measurement

During operation, the accuracy of oximetry readings can be affected by the following factors:

2.5.1 Instrument performance depends on the pulsatile character of the artery. The measurement would not be considered reliable and accurate if the following conditions are present during measurement:

- The patient is in cardiac arrest or shock
- Venous pulsations
- Low temperature of hand
- Have taken vascular activity medicine
- Anemia
- carboxyhemoglobin
- methemoglobin
- methylene blue
- Indocyanine green

2.5.2 Instrument performance depends on the wavelength absorption for oxyhemoglobin and deoxyhemoglobin. If there are substances absorbing the same wavelength, this would induce false or low SpO₂ values. The following may affect these values:

- carboxyhemoglobin
- methemoglobin
- methylene blue
- Indocyanine green

2.5.3 Extremely high illumination could affect the SpO₂ measurement. Use a semi-translucent or opaque cover to shield the sensor.

2.5.4 Other factors

- a) High-frequency electrosurgical interference from external devices, including defibrillators;
- b) Placement of a sensor on an extremity that currently has been installed a blood pressure cuff, arterial catheter, or intravascular line;
- c) The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia;
- d) An arterial occlusion proximal to the sensor.



Warnings

- Use only SpO₂ sensors provided by manufacturer. Other SpO₂ sensors may cause improper performance.
- Do not use an SpO₂ sensor with exposed optical components.
- Excessive patient movement may cause inaccurate measurements.
- The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia.
- Tissue damage can be caused by incorrect operation or misusing sensor; for example, by wrapping the sensor too tight. Inspect the sensor site to ensure the skin's integrity and the adhesion position of the sensor is correct. More frequent inspection should be taken if necessary.
- Loss of pulse signal can occur in any of the following situations:
 - a) The sensor is too tight;
 - b) There is excessive illumination from light sources such as a surgical lamp, a bilirubin lamp, or sunlight;
 - c) A blood pressure cuff is inflated on the same extremity as the one to which an SpO₂ sensor is attached.

CHAPTER 3 DETAILED OPERATION

3.1 Display

The handheld pulse oximeter uses an OLED display for a readout. It can display the SpO₂ and pulse rate (PR) value, as well as a pulse column and SpO₂ waveform.

There are three display modes shown in Fig 3.1. The first figure is pulse column display mode. The second figure is filled waveform mode. The third figure is line waveform mode indicating SpO₂% trend.

(1) Pulse column display mode. (2). Filled waveform mode. (3). Line waveform mode




Fig.3.1 Three display modes

- **SpO₂:** Percent oxygen saturation value displayed above is 98%
- **PR:** Pulse rate value displayed is 72 bpm
- **Pulse column:** This is used for signal identification and quality indication during motion and low signal to noise situations. The bar rises and falls with the pulse, its height indicating signal quality. When the bar is very low, the SpO₂ and pulse rate values may be suspect.
- **Signal strength:** Indicates arterial pulse signal strength and may be used as a diagnostic tool during low perfusion for the accurate prediction of illness severity. The bar is highest when the quality of the perfusion state is best and low when the perfusion is poor.
- **PR tone modulation:** Beeps in sync with the patient's pulse, even under most challenging patient motion conditions, except when in the (3) mode.

3.2 Key Functions

1) Power On

Press the Function key (left key) to power on.

The oximeter will power off automatically under the "information display mode" or "trouble mode" for 30 seconds.

2) Definition of key

There are two keys in the oximeter: "Function key (left key)" and "Setting key (right key)".

"Function Key": Located at the left of the oximeter, this acts as a Power On switch when the unit is in the Off condition. When the unit is On, it acts as a Function key.

"Setting key": Located at the right of the oximeter, it has no function when the power is off. When the unit is On, it acts as a Setting key.

3) Definition of Key Press

There are three ways to press the key:

Press: press the key quickly, the duration time should no more than 1 second.

Double press: Two-time continuous press, the time between the two press actions should be no more than 0.5 second.

Extended press: Press the key for more than 1.5 seconds.

3.3 Mode Introduction

Power On - Measurement and operation can be done normally. There are three display modes: "Measure mode", "Information display mode" and "Trouble display mode".

1) Display modes

"Measure mode": the sensor is plugged in correctly, and the finger is properly in the sensor, the oximeter is in the measurement mode for both SpO₂ and PR.

"Information display mode": The sensor is not plugged into the oximeter or the sensor is plugged in the oximeter but the finger is not in the sensor, "Probe off" or "Finger out" will be displayed respectively, and it will automatically power off if either information display lasts for more than 30 seconds.

"Trouble display mode": In the failure state, the oximeter will display error information, and will automatically power off if the information display lasts for more than 30 seconds. For error information details and definitions, please refer to chapter 4.2.

2) Display Mode setting

The Display Mode setting can only be enabled under "Measure mode" and "Information Display mode". Pressing the Function key (left key) or Setting key (right key), the display mode may be changed sequentially between "Filled Waveform", "Line Waveform" and "Pulse Column Waveform". Once you choose a certain mode, the oximeter will continually display in this mode until changed by the user.

3) Setting Mode

Enter the Setting Mode under the "Measure mode", you can only set Brightness and Patient ID. Enter the Setting Mode under the "Information Display mode", you can set Brightness, Patient ID, Date and Time. Current parameters and data will be displayed at the top right corner which is used to display the "PR". The various parameters display title is as follows:

Br – Brightness (range:1-7)
ID – Patient ID (range:1-10)
Y – Year (range:2000-2099)
M – Month (range:1-12)

D – Day (range:1-31)
H – Hour (range:0-23)
M – Minute, (range:0-59)
S – Second (range:0-59)

a) Enter “Setting mode”

Press the function key (left key) for more than one second (extended press), the oximeter will enter into “Setting mode” and you will find a parameter item which includes its title and value on the top right corner of the display (refer to Fig.2.5).

After entering the setting mode from “Information Display mode”, continually press “Function key (left key)”, the current parameter item will be changed sequentially in the following order: Br(brightness) -->ID(patient’s ID) -->Y(year) -->M(month) -->D(date) -->H(hour) -->M(minute) -->S(second) -->Br(brightness)-->.

Notice: You cannot set the date and time when in the measuring course.

b) Save and exit from “Setting mode”

- ① Under the “Setting mode”, press “function key (left key)” to select the desired parameter item.
- ② Then press “Setting key (right key)” to adjust the value.

Each time you press the “Setting key (right key)”, current parameter value will be added by 1 unit sequentially.

Double press the “Setting key (right key)” and the current parameter setting will be added by 9 units sequentially.

- ③ Press “function key (left key)” to select the next desired parameter item. Then redo step ② to adjust the value. You can continually redo step ① to ② until no parameter’s setting need to be changed.
- ④ To finish setup, **fleety press both left and right keys together to confirm**. Then the modifications under the “Setting mode” will be saved, at the same time, the system will exit from “Setting mode”.

Notice:

- To confirm the modification, the act of pressing both keys simultaneously must be fleet. Do not press them for long time.
- The date and time adjustment only can be done in the “Setting mode” under “Information Display mode”.

c) Cancel and exit from “Setting mode”

Under the “Setting mode”, double-press the Function key (left key), the modification under “Setting mode” will be cancelled and simultaneously, the system will exit from “Setting mode”.

If there is no operation under “Setting mode” for more than 30 seconds, the system will exit from “Setting mode” automatically.

3.4 Data replay and transmission

The oximeter can record SpO₂ and PR value for more than 72 hours, and can analyze records one by one. You can transfer the history data to a PC by using “MEDVIEW” software and an attached data cable. As for detailed setup and operation, please refer to the “MEDVIEW” operator’s manual.

Notice: The mode of its storage is cycle store, it only can save the latest 72-hour records.

3.5 Alarm

Alarm: Technical alarm and physiological alarm.

Technical alarm: finger out, probe off, power low and error code.

In the situation that the finger is not inserted correctly or the connection state of the probe is not good results in failure of measurement, "Finger out" or "Sensor off" may be displayed on the normal screen.

When battery power is lower than $2.3 \pm 0.1V$, the sign  will flicker in its display area. Replace the batteries as soon as possible.

In the failure state, the oximeter will display error codes, and will automatically power off if the error code display lasts for more than 30 seconds. For the details and definitions on error, please refer to chapter 4.2.

Physiological alarm: SpO₂ and PR

If the measured SpO₂ and/or PR value is beyond the default limit, alarm will be activated, and the corresponding value will flash with audible sound.

By default:

SpO₂: The upper limit: 100%

The lower limit: 90%

PR: The upper limit: 100bpm

The lower limit: 60bpm

Notice: During the alarm is issued, press the function button, you can silence the alarm for 30 seconds. When press it again, the screen switches between pulse column display mode, filled waveform mode and line waveform mode.

Alarm priority:

There are three-level priorities for selection.

High priority: indicates the patient is in the very dangerous situation.

Medium priority: indicates the warnings should be paid attention to.

Low priority: indicates the technical alarm caused by the device itself.

Alarms of the oximeter include technical and physiological alarms. All the three priorities are divided by built-in module and can not be changed by user.

Assignment of priority:

Alarm \ priority	Event	Display	Audible sound
High	The SpO ₂ is beyond the limits	The SpO ₂ value is flashing	"Di- Di – Di ---- Di - Di", "Di- Di – Di ---- Di - Di" once circularly every 3 seconds
Medium	The pulse rate is beyond the limits	The pulse rate value is flashing	"Di - Di - Di", once circularly every 5 seconds
Low	The probe or finger is not inserted	Finger out Probe off	"Di", once circularly every 20 seconds

CHAPTER 4 MAINTENANCE AND REPAIR

Warning: The advanced circuit inside the oximeter does not require periodic calibration and maintenance, except replacing the batteries.

Don't open the cover of oximeter or repair electronic circuits. Its open will cause the damage of the device and the annulment of the guarantee.

4.1 Maintenance

It is very important for user to perform daily maintenance of oximeter and parts in order to maintain its function and appearance. Disinfection procedures may be performed with the use of the below mentioned cleaner/disinfectants. Failure to perform these procedures may result in invalidating the warranty. Local disinfection protocols will apply.

The external surface of the oximeter can be cleaned by wiping with a damp cloth. Do not submerge the oximeter in any solution at any time. To do so will void the warranty

Use the following permitted solutions:

- Ammonia (diluted)
- Glutaraldehyde
- 10% Bleach solution
- Mild soapy water (diluted).
- Do not use the following cleaners:
 - Any kind of scrubbing or scouring solution
 - Acetone
 - Alcohol-based cleaners

Battery maintenance

Please take out battery before cleaning the oximeter.

Remove the batteries if you are not going to use the oximeter for a long time.

4.2 Troubleshooting

a) Error Definitions

Err 1: program memory damaged.

Err 2: data memory damaged.

Err 3: sensor Red Emission Diode damaged.

Err 4: sensor Infrared-ray Emission Diode damaged.

Err 5: sensor Infrared-ray Receipt Diode damaged.

Err 6: exterior crystal oscillator damaged.

Err 7: sensor emission diode or receipt diode damaged.

Err 9: real time clock damaged.

Err 10: EEPROM chip damaged.

b) Possible problem and corresponding resolution

Problems	Possible reason	Solution
SpO ₂ or PR cannot be displayed normally	<ol style="list-style-type: none"> 1. Finger is not plugged correctly 2. Patient's Oxyhemoglobin value is too low to be measured 	<ol style="list-style-type: none"> 1. Retry by plugging the finger 2. Attempt several times to obtain a reading, If you are sure that no problem exists, obtain further clinical examination
SpO ₂ or PR display is unstable	<ol style="list-style-type: none"> 1. Finger might not be plugged deep enough 2. Finger is trembling or patient is moving continually 	<ol style="list-style-type: none"> 1. Retry by plugging the finger 2. Urge the patient to remain still
The Oximeter can not be powered on	<ol style="list-style-type: none"> 1. Battery power may be inadequate or not installed 2. Batteries might be installed incorrectly 3. The Oximeter might be damaged 	<ol style="list-style-type: none"> 1. Please replace batteries 2. Please reinstall the batteries 3. Contact local customer Technical Service
"Probe off" displayed on screen	<ol style="list-style-type: none"> 1. The sensor is not connected 2. The connection between the Probe and Oximeter is loose 	<ol style="list-style-type: none"> 1. Connect the sensor 2. Please check if the probe was connected with oximeter correctly

4.3 Warranty and Repair**4.3.1 Service Method**

Service support: Our company will offer user telephone and e-mail support as well as in-house repair at our facility

Parts replacement: Our company will replace parts, accessories, free of charge during the warranty period.

Our company will update the system software free of charge.

4.3.2 Exempt and limitation:

- a) Our company isn't responsible for such damage caused by force majeure.
For example: fire, thunder flash, flood, cyclone, hail, earthquake, house collapse, commotion, plane failing and traffic accident, deliberate damage, lack of fuel or water, labor and capital bother, strike and stop-working etc.
- b) No-service offer
The corresponding charge and insurance charge of disassembling, refurbishing, repackaging and moving the oximeter or the part of it.
The damage caused by the third company not commended by our company to adjust, install, replace the parts of the oximeter.
The damage and failure caused by user or its representative doesn't comply with the operator's manual.
- c) The oximeter is installed or connected with such external device without our company permission as printer, computer, netline and lead to oximeter failure. Our company will charge for the maintenance.
- d) Responsibility limitation
During the period of maintenance contract validity, if user changes the parts manufactured by other manufacturers without our company permission, our company is entitled to stop contract.

4.3.3 User Guarantee

- a) Please read user manual carefully before operation.
- b) Please operate and make daily maintenance as request of manual and guarantee.
- c) Power supply and environment.

4.3.4 No-guarantee principle

- There is no-dispelled smut and not-original mark in the crust.
- There is physical damage on oximeter and its accessory.
- There are liquid leftover and eyewinker on oximeter, which may lead to short circuit and plug board failure.
- All the probe and accessories belong to consumption and beyond free change range.
- Such damage of probe caused by mechanical force doesn't belong to free change range.
- During measurement of SpO₂, principle leads to measuring value difficultly or inaccurate measurement.
- Not-original package leads to damaging oximeter during transportation
- Not our company professionals or authorized personnel disassemble oximeter and lead to oximeter failure.
- Not carefully read manual and so wrong operation lead to oximeter damage and failure.

4.3.5 User's Special Request for Guarantee Time

Our guarantee constitution for oximeter complies with electronic product after-sale service standard regulated by national laws. We regulate the guarantee time of hoist board is one year and all the accessories are three months. If users request the guarantee time beyond our regulated guarantee time, we should take it into consideration. Because electronic product has such character of quick changing, for such user asking more than three years guarantee time, our company will not sell oximeter parts during maintenance. Our company will upgrade oximeter or change new maintenance methods, for this, we charge the lowest price for new oximeter with user permission.

4.3.6 Repackage

- Take all the accessories and put them into plastic cover
- Try to use original package and packing material. User will be responsible for such damage caused by bad package during transportation.
- Please offer guarantee list and copy of invoice to standby with the period of guarantee.
- Please describe failure phenomenon in detail and altogether offer oximeter.

Storage and Transportation

Storage: Storage Temperature -20°C~55°C, Relative Humidity <93%, no condensation.

Transportation: Transport by airline, train or vessel after packing according to request.

Package

Pack the product with the hard bag, and put the foam between the inner box and the carton to alleviate the shake.

APPENDIX A Specifications

Notices:

- Specifications may be changed without prior notice.
- The circuit diagrams, the list of components, the illustration of diagrams, and the detailed rules of calibration, are provided exclusively to professional personnel authorized by our company.

SpO₂

Display Range: 0%-100%

Measurement Range: 70%-100%

Resolution: 1%

Accuracy: 80%-100% : $\pm 2\%$; 70~79% $\pm 3\%$; <70% : unspecified

Data update time: < 15 s

Probe LED Specifications:

	Wavelength	Radiant Power
RED	660 \pm 2nm	1.8mW
IR	940 \pm 10nm	2.0mW

Pulse(Heart) Rate

Display Range: 0-254 bpm

Measurement Range: 30-235 bpm

Resolution: 1bpm

Accuracy: ± 2 bpm or $\pm 2\%$

Alarm

By default: SpO₂: The upper limit: 100%

The lower limit: 90%

PR: The upper limit: 100bpm

The lower limit: 60bpm

Other Alarms: Probe off, Finger out, Low battery

Modes: Visual and audible

Display

Type: OLED, double color

Parameters: SpO₂, PR, Pleth waveform, Pleth bar

Mode: 3 display modes.

Record

Patient ID: 10 patients

Record time : Up to 72 hours (cycle store)

Data transmission

Transmission method: Cable Transmission

Environmental

Operating Temperature: 5°C~40°C

Storage Temperature: -20°C~55°C

Operating Humidity: ≤80%RH, no condensation

Storage Humidity: ≤93%RH, no condensation

Classification per IEC60601-1

Type of protection: Internally powered equipment

Degree of protection: Type-BF applied Par

Mode of Operation: Continuous

Safety: IEC Standard 60601-1

Dimensions: 130 x 44 x 20 mm (Length×Width×Height, not contain the probe)

Weight: 120g (with alkaline batteries)

Power supply

Type: 2 AAA Alkaline Batteries

Operation time: About 16 hours of typical operation

APPENDIX B The declare of manufacturer

Guidance and manufacture's declaration - electromagnetic immunity for Equipment and Systems that are not Life-Supporting

Guidance and manufacture's declaration - electromagnetic immunity

The Handheld Pulse Oximeter is intended for use in an electromagnetic environment specified below. The customer or the user of the Handheld Pulse Oximeter should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance Level	Electromagnetic environment guidance
Electrostatic Discharge (ESD) IEC61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floor are converted with Synthetic material, the relative humidity should be at least 30%
Electrical fast transient/ burst IEC 61000-4-4	Not applicable		
Surge IEC 61000-4-5	Not applicable		
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable		
Power frequency magnetic field IEC 61000-4-8	3 A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

Guidance and manufacture's declaration - Electromagnetic Immunity for Equipment and Systems that are not Life-Supporting

Immunity Test	IEC 60601 Test level	Compliance Level	Electromagnetic environment guidance
Conducted RF IEC 61000-4-6	Not applicable		<p>Portable and mobile RF communications equipment should be used no closer to any part of the handheld pulse oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ <p>Where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacture and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range ^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol.</p> <div style="text-align: right;">  </div>
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5 GHz	3V/m	

NOTE1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE2 These guideline may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base situation for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If abnormal performance is observed, the additional measures may be necessary, such as reorienting or relocating the handheld Pulse Oximeter.

^b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Handheld Pulse Oximeter

The Handheld Pulse Oximeter is intended for use in an electromagnetic environment in which radiated RF disturbance are controlled. The customer or the user of the Handheld Pulse Oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the Handheld Pulse Oximeter as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150KHz to 80 MHz $d = \left[\frac{3.5}{V_i} \right] \sqrt{P}$	80MHz to 800 MHz $d = \left[\frac{3.5}{E_i} \right] \sqrt{P}$	800MHz to 2.5 GHz $d = \left[\frac{7}{E_i} \right] \sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.689	3.689	7.379
100	11.667	11.667	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic interference is affected by absorption and reflection from structures, objects and people.

The Manufacturer Information

Beijing Choice Electronic Technology Co., Ltd.

Add: Bailangyuan Building B, Rm. 1127-1128, Fuxing Rd., A36, Beijing, China

Tel: 0086-10-88204188

Fax: 0086-10-88203199

E-mail: chenlei@choicemmed.com

The ER Representation Information

Shanghai International Trading Corp. GmbH (Hamburg)

Add: Eiffestrasse 80, 20537 Hamburg Germany

Tel: 0049-40-2513175

Fax: 0049-40-255726

E-mail: antonsissi@hotmail.com

